



The SKYTRON 6001 Surgical Table is an electro-hydraulic surgical table with a remote control and C-arm positioning to meet all imaging requirements. The stability of the table and the rigidity of the top provide assurance that the patient position will be maintained and desired.

Features

- C-Arm Compatible Table Top
- Wired Hand Control
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- 15-foot power cord with hospital grade plug
- Standardized side rails
- Four-point, push-button operated braking system for safety and stability
- 30° lateral tilt
- Trendelenburg and reverse capabilities of 30°
- Weight Capacity: 500 lbs.



Specifications

General Specifications

Power Requirement: 120 VAC, 60 Hz.

Weight Capacity: 500 lbs.

Height Minimum: 28".

Height Maximum: 45".

Width: 20".

Length: 75".

Lateral Tilt: 30°.

Trendelenburg: 30°.

Top Rotation: 180°.

