



Uncommon flexibility makes the Skytron 6500 HD general surgical table unequaled for general-purpose surgery. The exclusive 180° rotating top feature of the Skytron 6500 general surgical table allows unobstructed C-arm and microscope coverage. Additionally, it eliminates conflicts between surgeon and anesthesiologist during head-end procedures and provides ample legroom for surgeons who prefer to sit during procedures.

## Features

- Lift and fully articulate patients weighing up to 1000 lbs.
- Length with Headrest is 79 inches.
- The tabletop can rotate 180 degrees.
- Built-in Powered Kidney Bridge.
- Removable powered leg section.
- Self-leveling, four-point brake system.



# Specifications

## Dimensions

Length: 65.5" (166.5 cm)  
Length (with head rest): 79" (201 cm)  
Width: 20" (50.4 cm)  
Width (with side rails): 21.75" (54.8 cm)  
Maximum Permissible Load (lift): 1000 lbs (454 kg)  
Maximum Permissible Load (Articulation): 1000 Lbs (454 kg)

## Hand Control Adjustments

Height: 26" - 47" (65.5 - 118.4 cm)  
Lateral Tilt (left/right): 30°  
Trendelenburg/Reverse Trend:  
30°/30°  
Power Kidney Bridge: 5" (12.6 cm)  
Head Section (up/down) 60°/90°  
Back (up/down): 90°/40°  
Leg (up/down): 20°/150°  
Top Rotation: 180°  
Flex/reflex: 40°/90°  
Beach Chair: Available

## Electrical

Power Requirements: 120 VAC, 60 HZ 300 watts  
Current Leakage: Less than 100 microamps  
Power Cord: 15 Feet with standard duplex connector