



The **Skytron 6700B** is a general surgical table that provides 1,200 lb. lift and 1,000 lb. articulation capacity with removable back and leg sections. These removable back and leg sections can optionally be carbon fiber to further extend the imaging window. The 6700B surgical table features an “ultra-low” position of only 23-inches to help combat the challenges of patient transfer.

Features

- 210° rotating top.
- Auto beach chair position.
- 30° tilt, trendelenburg, reverse trendelenburg.
- 4-point auto-leveling braking system.



Specifications

Dimensions

Height: 23 - 41" (58.4 - 104.1 cm)
Width: 19.75" (50.1 cm)
Width (With side rails): 22" (55.9 cm)
Length: 65.5" (166.3 cm)
Length (With headrest): 79" (200.7 cm)
Weight: 900 lbs (408.2 kg)

Movements

Height: 23 - 41" (58.4 - 104.1 cm)
Lateral tilt (left/right): 30°
Trendelenburg: 30°
Reverse Trendelenburg: 30°
Power Kidney Bridge: 6" (15.2 cm)
Back up/down: 90°/40°
Leg up/down: 20°/100°
Flex/Reflex: 40°/90°
Top Rotation: 40°

Max. Load

Weight Limit (Lift): 1,200 lbs. (544.3 kg)
Weight Limit (Articulate): 1,000 lbs. (453.6 kg)

Battery

Optional

Return to level

Yes

