



The **Skytron 6701 Hercules** surgical table has a 1,200 lb. lift and 1,000 lb. articulation capacity, removable back and leg section, 210° top rotation, 23" low table height, auto beach chair position, and more! The 210° table top rotation provides great C-arm flexibility and helps eliminate conflicts between the surgeon and anesthesiologist during head-end procedures.

Features

- 1,200 lb. lift and 1,000 lb. articulation capacity.
- 30° patient articulation including tilt, trendelenburg, and reverse trendelenburg.
- 23" low height adjustment.
- Full body imaging coverage.
- Automatic 4-point, self-leveling brake safety system.
- Optional bariatric table accessories.



Specifications

Dimensions

Height: 23-41" (58.4-104 cm)
Width: 19.75" (50 cm)
Width with side rails: 22" (56 cm)
Length w/o head rest: 65.5" (166.5 cm)
Length with head rest: 79" (201 cm)
Weight: 900 lbs (410 kg)

Movements

Height: 23-41" (58.4-104 cm)
Lateral Tilt (left/right): 30°
Trendelenburg: 30°
Reverse Trendelenburg: 30°
Power Kidney Bridge: 6" (15 cm) Lift
Back up/down: 90°/40°
Leg up/down: 20°/100°
Return to level (anatomically correct return to horizontal position from articulated positions)
Top Rotation: 210°
Flex/Reflex: 40°/90°
Base braking (lock/unlock)
Beach Chair

Max. Permissible Load

1200 lb. (545 kg) Lift
1000 lb. (455 kg) Articulation

