



The 3000 Imaging Table with 83" Extended Carbon Fiber Top and 500 lb. weight capacity is ideally suited for pain management, etc. The 3000 Imaging Table with 83" Extended Carbon Fiber Top and 500 lb. weight capacity is ideally suited for pain management, vascular and angiography procedures, while providing an unobstructed imaging window of 27.7" X 65".

Features

- 500 lb. Capacity
- 83" Extended carbon fiber top
- 65" imaging area
- 20 degree lateral tilt, Trendelenburg, Reverse Trendelenburg
- Height adjustment 28" to 39"
- Narrow base design for C-arm access
- One touch pendant control operated at any point around periphery of table.



Specifications

Dimensions

Length: 83" (2,100 mm)

Width: 19 5/8" (500 mm)

Adjustments Via Hand Control

Height: 28"-39" 710-990.6 mm

Lateral tilt (left/right): 20°

Trendelenburg: 20°

Reverse Trendelenburg: 20°

Max Permissible Load

500 lb (227kg) Lift

500 lb (227kg) Articulation

