



The **Skytron Elite 6500** series surgical tables are electro-hydraulically operated, general purpose surgical tables. The electro-hydraulic positioning functions are operated by the hand-help, push button, pendant controller. On the controller, you have the options for Trendelenburg, Lateral tilt, Back section, Elevation, Leg section, Flex/reflex, Kidney life, Return to level, and the floor lock/brake system. There are manual controls that provide positioning for the head section, table top rotation, emergency brake release, and leg section removal.

Features

- 500 pound weight capacity
- 180° rotating top
- Detachable head section
- Provides full body imaging
- Easy-to-use pendant control



Specifications

Dimensions
Height: 26 - 47" (66 - 119.4cm)
Width (Without side rails): 20" (50.8 cm)
Width (With side rails): 21.75 (55.24 cm)
Length: 76" (193.04 cm)
Patient Weight Capacity: 500 lbs (226.8 kg)

Trendelenberg / Reverse Trendelenberg
Trendelenberg: 30°
Reverse Trendelenberg: 30°

Head Section
Raise: 0-15°
Lower: 0-90°

Tilt
Lateral Tilt: ±30°

Back Section
Down: 40°
Up: 90°

Seat Section
Flex: 0-70°
Reflex: 0-90°

Leg Section
Down: 105°
Up: 20°

Kidney Bridge
Elevation: 5" (12.7cm)

Rotation
 210°

