



The **Skytron Essentia 1700** is engineered to deliver quality, value and long term durability. The Skytron 1700 is a general purpose surgical table is built to the same rigorous standards as Skytron specialty surgical tables, equipped to provide full body imaging windows and 500 lb. capacity and menu of essential table features. The 1700 is equipped with a removable 15 ft. power cord, 4 conductive casters, 4 point braking system with push button operation, including manual override.

## Features

- 500 lb. Capacity.
- Full Body Imaging.
- 180° Top Rotation.
- One-Touch Pendant Control.
- Line Power.
- Manual Emergency Brake Release.



## Specifications

### Dimensions

**Height:** 29 to 39" (73.7 to 99 cm)

**Width:** 22" (56 cm)

**Length (w/o head rest):** 65.5" (166.4 cm)

**Length (with head rest):** 76" (193 cm)

**Weight:** 530 lbs (240 kg)

### Table Movements

**Height:** 29 to 39" (73.7 to 99 cm)

**Lateral Tilt (Left/Right):** 25°

**Trendelenburg:** 25°

**Reverse Trendelenburg:** 25°

**Back Up/Down:** 60°/30°

**Leg Down:** 90°

### Weight Capacity

**Capacity:** 500 lbs. (227 kg)

