



The **SmartVest SV2100** is a high-frequency chest wall oscillation (HFCWO) device. The SmartVest airway clearance system is designed to deliver high-frequency chest oscillation to promote airway clearance and improve bronchial drainage. HFCWO has been proven to be an effective therapy for a wide range of conditions. And when patients use the simple, comfortable SmartVest HFCWO system as intended, the changes can be life-transforming. The portable air-pulse generator offers unique and patented features expressly designed to enhance ergonomic comfort and lifestyle convenience, helping patients get better by making it easier and more efficient to do therapy consistently. The air pulse generator on the SmartVest creates air pulses that are transmitted to the chest via the connection hose and vest. The generator is designed to operate in either the upright or horizontal position.

Features

- Washable SmartVest garment. (sold separately)
- Removable air-bladder system.
- Programmable, multi-position, ultra-quiet SmartVest SV2100 generator.
- Connection hose Air/Oxygen Mixer.
- Available in different sizes.

