



The **Sonosite HFL 50x** is a 15-6 MHz multi-frequency, broadband, 50mm linear array transducer for small parts, breast, musculoskeletal, and nerve exams. Supports 2D, SonoMBe (Proprietary software that enhances the system's ability to aid needle visualization), M-mode, Color Power Doppler, Velocity Color Flow Doppler and Pulsed Wave Doppler. Maximum depth is 6 cm.

## Features

- 15-6MHz
- Linear Array Transducer
- Scan Depth: 6 cm
- Small Parts applications

