



With the ability to produce hard copy images in just 10 seconds, Sony's UP-960 Black and White Video Graphic Printer is ideal for medical black and white applications such as ultrasound, mobile c-arm, and cardiac cathode lab. Features: Highlights True 256 gray scale gradations. Video loop through. UL2601-1 listed for medical use. Worldwide power and video signal capability. 10 Second print time

Features

- True 256 gray scale gradations
- Video loop through
- UL2601-1 listed for medical use
- 10 Second print time
- Worldwide Power Supply
- Small Footprint



Specifications

Dimensions	W x H x D 316mm x 132.5mm x 305mm
Printer Gradation	256 gray scale gradation
Power Requirments	100V to 240V, 50/60 Hz
Print Speed	10 seconds / 150mm NTSC
Printer Driver Software	No
Printing Method	Thermal graphic printing
Printer Resolution	No
Saftey Compliance	No
Warranty	1 Year

