



The Axia Scope Cabinets can be ordered in several different models depending on the amount of scope hangers your facility needs. These scope cabinets are a fully welded and stainless steel construction with glass hinged doors with louvered vents. There is a lower rack in the middle of the cabinet to keep the scopes in alignment. The standard features of the Axia scope cabinet are listed below

Features

- Glass doors with louvered vents.
- Scope hangers ranging from 9-22 hangers.
- Fully welded stainless steel construction.
- Lower middle alignment rack in cabinet.
- Removable drain pan.
- Padded back wall to prevent damage to the scopes.



Axia Scope Cabinet Specifications

Models and Dimensions

- AX-SC9: (HxWxD) 92" x 36" x 18" with 9 Scope hangers.
- AX-SC22: (HxWxD) 92" x 42" x 24" with 22 Scope hangers
- AX-SC9PV: (HxWxD) 92" x 36" x 24" with 9 Scope hangers and a power vent.
- AX-SC22PV: (HxWxD) 92" x 42" x 24" with 22 Scope hangers and a power vent.

