



The **Steris Amsco 2532** is a single chamber washer and disinfector. It's designed to support smaller CSSD's and ambulatory surgery centers. The Steris 2532 utilizes a full coverage spray pattern with cross pattern spray nozzles to ensure complete coverage. This washer features a drop-down loading door that allows loading and unloading trays without the need of a transfer cart and an intuitive LCD touchscreen control panel which provides clear visibility and access to critical information.

Features

- Manual drop-down doors with interior and exterior glass to allow the operator to view the chamber.
- Stainless-steel pump that provides circulation using a 3 HP motor.
- Drains between rinses to minimize residue remaining on instruments.
- Removable stainless-steel filter
- 3 Chemical injection pumps
- High-capacity air blower that delivers 240°F (116°C) of fresh air through unit piping. The fresh air is drawn through a HEPA filter.
- 5.7" (11.5cm) color touchscreen



Specifications

Exterior Dimensions

Height: 78" (198.1 cm)
Width: 28.5" (72.4 cm)
Depth: 28.5" (72.4 cm)
Weight: 540 lbs (245 kg)

Interior Dimensions

Height: 25.5" (64.8 cm)
Width: 25" (63.5 cm)
Depth: 25.25" (64.2 cm)

Utility Requirements

Hot Water: 1/2" NPT
Cold Water: 1/2" NPT
Ventilation: 4" (10.2 cm) O.D.
Drain: Recommended minimum 3" (7.6 cm) drain outlet
Pure Water: 1/2" NPT

Power

208 V, 60 Hz, 3-Phase, 60 Amp
380-400 V, 60 Hz, 3-Phase, 30 Amp
460-480 V, 60 Hz, 3-Phase, 30 Amp
380-415 V, 50 Hz, 3-Phase, 30 Amp

Display

Type: LCD
Size: 5.7" (11.5 cm)

Load Height

29-3/4" (75.6 cm) above finished floor

Chamber and Cycle Times

Cycle Time: 30 Minutes
Rack Capacity: 6
Cycles Per 8 Hour Day: 16
Trays Per 8 Hour Day: 96
Total Water Used Per Cycle: 21.2 Gallons (80.25 Liters)

