



The AMSCO 5052 Single-Chamber Washer/Disinfecter is designed for use in the cleaning and intermediate-level disinfection of soiled reusable utensils, trays, bedpans and urinals. It can also be used for rubber and plastic goods, simple hard-surfaced rigid surgical instruments (such as forceps and clamps) and other similar and related articles found in healthcare facilities.

Features



- Color touch screen microprocessor control system
- Audible warning system
- Space saving, power vertical sliding doors
- Eight factory-loaded cycles. Each cycle effective, fast and energy efficient. Each cycle may be customized (within allowed range of parameters) to meet specific operating requirements
- Chamber illumination by long-life, energy efficient LED light
- Compact space saving footprint
- Conforms to seismic building code requirements

Specifications

External Dimensions

Height: 80 in (2032 mm)
Width: 42 in (1067 mm)
Length: 31 in (787 mm)

Interior Chamber Dimensions

Height: 26.25 in (667 mm)
Width: 26.5 in (673 mm)
Length: 26 in (660 mm)

Load Height

31 in (787 mm) above the finished door

