



The **Steris Amsco 3085 R** surgical table is available with ACT ENABLED® electronic interface circuits that work seamlessly with other technology in your state-of-the-art OR. This integration system agnostic technology enables seamless communication among virtually any vendor's equipment.

Features

- Automatic conversion of functions for reverse patient positioning, regardless of headrest location means commands are automatically interpreted according to the selected orientation.
- Bieri hydraulics provide ultra-precise tabletop articulation for specialties such as urology and neurology. Tabletop movement is smooth and accurate through the full range of positions.
- With the addition of the orthopedic extension, the Amsco 3085 R surgical table pulls traction and provides positioning flexibility for hip and lower extremity orthopedic procedures.
- T-shaped pedestal design and reverse orientation provide optimal access to the patient by the entire surgical team, regardless of the procedure.



Specifications

Dimensions

Table Length: 76" (193 cm)

Table Width: 20" (51 cm)

Table Height Range: 27" to 44" (68.6-111.8 cm)

Patient Weight Capacity:

- 1,000 lbs (454 kg) - normal orientation
- 500 lbs (227 kg) - side tilt & reverse orientation

Hand Control Features

Normal/reverse orientation button

Return to level button

Table Movements

Trendelenburg / Reverse: 25° / 25°

Lateral Tilt (left/right): 18° / 18°

Head Section: +90° to -90°

Back Section:

- +55° to -25° (normal)
- +80° to -105° (reverse)

Leg Section:

- +80° to -105° (normal)
- +55° to -25° (reverse)

Upper Body Imaging Area

Normal orientation (28" x 16"): 448 sq in (2890 cm²)

Normal orientation with IA extender (36" x 16"): 576 sq in (3716 cm²)

Reverse orientation, head section reversed (45" x 16"): 720 sq in (4645 cm²)

Reverse orientation 36" carbon fiber fluoro extension (36" x 20"): 720 Sq in (4645 cm²)

Lower Body Imaging Area

Normal orientation (33" x 14"): 462 sq. in. (2980 cm²)

Normal orientation, head section reversed (45" x 16"): 720 sq. in. (4645 cm²)

Normal orientation, 36" carbon fiber fluoro extension (36" x 20"): 720 sq in (4645 cm²)

Normal orientation with Fem/Pop board (55" x 16"): 880 sq. in. (5678 cm²)

