



The AMSCO 3085 SP surgical table is available with ACT ENABLED® electronic interface circuits that work seamlessly with other technology in your state-of-the-art OR. This integration system agnostic technology enables seamless communication among virtually any vendor's equipment.

## Features

- Automatic conversion of functions for reverse patient positioning, regardless of headrest location means commands are automatically interpreted according to the selected orientation.
- Bieri hydraulics provide ultra-precise tabletop articulation for specialties such as urology and neurology. Tabletop movement is smooth and accurate through the full range of positions.
- With the addition of the orthopedic extension, the AMSCO 3085 SP surgical table pulls traction and provides positioning flexibility for hip and lower extremity orthopedic procedures.
- T-shaped pedestal design and reverse orientation provide optimal access to the patient by the entire surgical team, regardless of the procedure.



# Specifications

## Dimensions

**Table Length:** 76" (193 cm)

**Table Width:** 20" (51 cm)

**Table Height Range:** 27" to 44" (68.6-111.8 cm)

**Patient Weight Capacity:**

- 1,000 lbs (454 kg) - normal orientation
- 500 lbs (227 kg) - side tilt & reverse orientation

## Hand Control Features

**Normal/reverse orientation button**

**Return to level button**

## Table Movements

**Trendelenburg / Reverse:** 25° / 25°

**Lateral Tilt (left/right):** 18° / 18°

**Head Section:** +90° to -90°

**Back Section:**

- +55° to -25° (normal)
- +80° to -105° (reverse)

**Leg Section:**

- +80° to -105° (normal)
- +55° to -25° (reverse)

## Upper Body Imaging Area

**Normal orientation (28" x 16"):** 448 sq in (2890 cm<sup>2</sup>)

**Normal orientation with IA extender (36" x 16"):** 576 sq in (3716 cm<sup>2</sup>)

**Reverse orientation, head section reversed (45" x 16"):** 720 sq in (4645 cm<sup>2</sup>)

**Reverse orientation 36" carbon fiber fluoro extension (36" x 20"):** 720 Sq in (4645 cm<sup>2</sup>)

## Lower Body Imaging Area

**Normal orientation (33" x 14"):** 462 sq. in. (2980 cm<sup>2</sup>)

**Normal orientation, head section reversed (45" x 16"):** 720 sq. in. (4645 cm<sup>2</sup>)

**Normal orientation, 36" carbon fiber fluoro extension (36" x 20"):** 720 sq in (4645 cm<sup>2</sup>)

**Normal orientation with Fem/Pop board (55" x 16"):** 880 sq. in. (5678 cm<sup>2</sup>)

