



The **Steris Amsco 600** medium steam sterilizer meets the critical demands of maximum throughput with a minimum footprint. The Steris Amsco 600 is designed for pre-vacuum and gravity sterilization of heat and moisture-stable medical devices and their accessories in healthcare facilities. This autoclave safely sterilizes up to 375 lbs per cycle for efficient high-volume sterilization.

## Features

- Easy loading with automated vertical door.
- No side or rear service clearance required.
- Maximizes space with maximum throughput.
- Active pneumatic valves for long-life reliability.
- Integrated large-capacity, high-quality vacuum pump.
- Fully jacketed chamber to maximize heat distribution and reduce condensation.
- Dual Pressure and temperature sensors.



## Specifications

### Dimensions

**Height:** 89.1" (226.3 cm)  
**Width:** 39.4" (100.1 cm)  
**Depth:** 51.7" (131.3 cm) or 63.4" (161 cm) or 75.3" (191.3 cm)  
**Floor Space:** 14.1 ft<sup>2</sup> (1.3 m<sup>2</sup>) or 17.3 ft<sup>2</sup> (1.6 m<sup>2</sup>) or 20.6 ft<sup>2</sup> (1.9 m<sup>2</sup>)

### Capacity

**Quantity of Wrapped Trays per Cycle:** 9, 12, or 15  
**Weight Per Tray:** 25 lbs (11.3 kg)  
**Total Weight Per Cycle:** 225 lbs (102.1 kg), 300 lbs (136.1 kg), or 375 lbs (170.1 kg)

### Thermal Printer

**Type:** High-resolution 8 dots/mm  
**Speed:** 25 lines/minute  
**Paper Size:** 2.25" (57 mm) wide thermal paper  
**Resistances:** Resistant to exposure from steam, alcohol, UV and visible light, oil, heat and water

### Environmental

**Pollution Degree:** 2  
**Installation Category (Overvoltage Category):** II  
**A-Weighted Sound Pressure Level:** 64.24 dB  
**IP (Ingress Protection):** Rating = IP20  
**Temperature:** 50 - 104°F (10 - 40°C)  
**Humidity:** 60% non-condensing, maximum  
**Altitude Range:** 0 - 6,500 ft (0 - 1981 m)

### Cycles

**Pre-programmed Cycles:** 6 pre-programmed cycles  
**Pre-programmed Test Cycles:** 3 pre-programmed test cycles

