



The AMSCO 7053L Single-Chamber Washer/Disinfector meets the demands for the highest-volume sterile processing departments. It is designed for cleaning and intermediate-level disinfection of soiled reusable utensils, trays, bedpans and urinals. It can also be used for rubber and plastic goods, simple hard-surfaced rigid surgical instruments (such as forceps and clamps) and other similar and related articles found in healthcare facilities.

The AMSCO 7053L Washer/Disinfector powers away soils with STERIS's exclusive Category 6 spray arm technology. Washer/disinfector is compliant with ANSI/AAMI ST15883-2 and ISO 15883-2 standards for thermal disinfection. Optional automation system minimizes handling time. USB port enables simple method for data management.

Features

- Productivity – 27 minute instrument cycle processes up to 15 DIN trays simultaneously
- Safe Operation – Easy to load racks minimize hands-on time with shelves that load from any side. Cool touch handles are staff-friendly
- Efficiency – Environmentally friendly design uses only 39 gallons of water per cycle



Specifications

Exterior Dimensions	(W x H x D) 42 x 80-3/4 x 38" (1067 x 2051 x 965 mm)
Chamber Dimensions	(W x H x D) 26-1/2 x 26-1/4 x 31.5" (673 x 667 x 800 mm)
Load height	(W x H x D) 31" (787 mm) above finished floor
Productivity	<p>A 21% increase in chamber size allows you to reprocess up to 50% more trays in an 8 hour day compared to the Vision® Single Chamber Washer/Disinfector</p> <p>Allows for flexibility when reprocessing a wide range of surgical instruments with multiple rack configurations</p>
Spray	<p>Full Coverage Spray Pattern combines our efficient cross pattern spray arm nozzles, impingement options and dual speed pumps.</p> <p>Innovative corner rack spray nozzles on each level, ensures no area of the tray is missed</p> <p>Providing maximum coverage, longer spray arms with nozzles at a 45 degree angle ensure no area is left untouched</p>
Drying	Reduces drying time with our unique four-sided air inflow while recycling the air for maximum efficiency
Trays	Holds up to 15 standard size trays
Tray Size	(W x L x D) 10" x 21.25" x 3.75
Easy of Use	<p>Cool touch handles allow for safe and immediate handling of loads following completion of cycle</p> <p>All racks are conveniently designed for easy load and unload from all sides</p>
Connectivity	Access critical cycle parameters at any time using the USB port for downloading
Control Panel	Increase staff confidence and safety with our intuitive LCD large touch-screen control panel which provides clear visibility and access to critical information