



The **Steris Amsco Dual Compartment Blanket Warmer** was designed to handle a large volume of patient warming needs. The blanket warmer has two different compartments one at 13 in height and the other at 36 in height. The warming cabinet is intended to be used to warm IV fluids, irrigation solutions, and linens and blankets. Each warming chamber is independently controlled to help fit your needs with efficient heating electronic door locks and electronic data reporting.

Features

- LED display controls for instant temperature reading
- adjustable shelves
- Door locks for security
- Temperatures can range from 90 - 160 degrees F Can warm and maintain temperatures for IV fluids, irrigation solutions, blankets and linens



Specifications

Dimensions

Outside dimensions: at 18" Depth: 74 3/4 x 30 x 20 1/2" (1899 x 726 x 521 mm)
Outside dimensions: at 24" Depth: 74 3/4 x 30 x 26 1/2" (1899 x 762 x 673 mm)
Inside Dimensions: at 18" Depth: Upper chamber: 13 1/2 x 24 x 16 7/8" (343 x 610 x 429mm)
Lower chamber: 36 1/2 x 24 x 16 7/8" (927 x 610 x 429mm)
Inside Dimensions: at 24" Depth: Upper chamber: 13 1/2 x 24 x 22 7/8" (343 x 610 x 581mm)
Lower chamber: 36 1/2 x 24 x 22 7/8" (927 x 610 x 581mm)

Technical Specifications

Door Option: Steel or dual-pane tempered glass
Compartments: 2
Total Capacity: at 18" depth = 11.6 cu. ft. (3.1 cu. ft. upper/8.5 cu. ft. lower)
Total Capacity: at 24" depth = 15.8 cu. ft. (4.2 cu. ft. upper/11.6 cu. ft. lower)
Installation: Free-standing or recessed
Mobile Base: Optional on 24" depth model only
Temperature Range: 90°F (32°C) to 160°F (71°C)
Electronic Door Locks: Optional