



The **Steris Amsco M70WC-EL** is a warming cabinet that is used for primarily heating surgical solutions in flasks with vented closures, but can also be used as a blanket warmer. The Steris Amsco M70WC-EL comes in two different models. One model is 18" deep and the other is 24" deep. Both of these units can be installed as open mounted or recessed into a wall. They also can either be electric or steam.

Features

- 2 Compartment, 18" or 24" depth.
- Steam or electric.
- Stainless Steel construction.
- Analog knob for heat adjustment.
- Adjustable stainless steel shelves in bottom compartment.
- Quiet operation.



Specifications

Dimensions (18" Depth)

Height: 74.75" (189.9 cm)

Width: 30" (76.2 cm)

Depth: 20.38" (51.7 cm)

Dimensions (24" Depth)

Height: 74.75" (189.9 cm)

Width: 30" (76.2 cm)

Depth: 26.38" (67 cm)

Heating System

Electric: Operates on 120-volt, 60-Hz power.

Steam: Has a thermostatically controlled heater valved and triggered for operation on steam from building supply delivered at 35-45 psig.

Heating Chamber

Upper compartment is insulated with 1" thick nominal fiberglass blanket. A fan circulates air within the chamber to provide even heat distribution. The door is a welded double wall construction with 1.5" thick insulation.

Lower compartment has adjustable stainless-steel shelves adjustable in 1" increments. Door of lower compartment is double wall construction with sound deadening between the walls.

Temperature Range: 100-180°F (38-82°C)

