



Steris Amsco Tabletop Digital Warming Cabinet: Model QDJ-03. The Steris Amsco QDJ-03 series of tabletop blanket warmers is a perfect solution to save space in your surgery center or department. The QDJ-03 is a cost-saving solution for blanket warmers in your facility versus the cost of a comparable new unit.

Features:

- Easy grip door handle.
- Temperature set dial, °f
- LED temperature window.
- Overtemp warning light.
- Rubber gasket in door closure.
- Chamber size is 23.5 x 22 x 13.5
- 21 x 23 in. Shelf Included.
- 4.5 foot heavy duty power cord.
- Power on/off rocker switch

Specifications:

Power	120V, 50.60 Hz.
Width	29.75"
Depth	28"
Height	25"
Weight	150 lbs.

