



The Reliance Synergy Washer/Disinfector is a mechanical washer equipped with an Eagle 3000 Stage 3 control system. The washer/disinfector is designed with nine adjustable cycles. One additional cycle is available for customized programming to meet specific operating requirements. Three injection pumps are provided with a standard washer/disinfector. One enzyme pump, one detergent pump and one lubricant pump (for Thermal Rinse phase) increase the flexibility of the customized cycles. Up to five injection pumps can be installed. Each preprogrammed cycle is equipped with PreWash, Enzyme, Wash, Rinse and Thermal Rinse phases.

Features

- 10 DIN trays Single-Chamber Double-Door Washer/Disinfector compliant with EN / ISO 15883
- Very good rack-to-chamber ratio of 96% (the closer to 100% the better the fit of the rack for the chamber)
- Steam or electric heating



Specifications

Size (W x H x L)

Exterior dimensions: 42 x 80 x 36-1/2" (1067 x 2032 x 927 mm)

Interior chamber dimensions: 26-1/2 x 26-1/4 x 26" (673 x 667 x 660 mm)

Load Height: 31" (787 mm) above finished floor

Dw rinse (CYCLE DURATION, MIN)

0-10

Rinse (FLUID RESERVOIR)

0-15

Rinse/lubricate (CYCLE DURATION, MIN)

0-10

Detergent wash (CYCLE DURATION, MIN)

0-15

Cold prewash (CYCLE DURATION, MIN)

0-15

Dry (CYCLE DURATION, MIN)

0-30

PUMP INTENSITY,HP (SIGNALS)

Dual-speed 7.5, 2

Max flow, gal/min (SIGNALS)

340.7 (90); 916 (242)

CONTROLLER (POWER NEEDED)

Microprocessor

HEAT SOURCE (CONTROL CONSOLE)

Steam or electric

WEIGHT, kg (lb) (DISPLAY)

544.3 (1,200.2)

