




The SurgiGraphic surgical table enables radiology-quality imaging for minimally-invasive vascular, orthopedic, or other procedures requiring fluoroscopy; lateral and longitudinal float, Trendelenburg and reverse Trendelenburg, lateral tilt, height adjustment, and the most useful combinations of these movements.

Features

- Unobstructed imaging from carotid to distal tibia - a full 83" - meaning a 360 degree view
- Higher quality imaging with less radiation energy due to minimized attenuation, scatter and edge effects
- Surgeon has complete control of every table movement with the Surgeon Palm Control
- Smooth longitudinal and lateral float



Specifications

General

Equipment Number TA-01

Unobstructed Imaging Length of Table 71" (180 cm)

Unobstructed Imaging Length with Carbon Fiber Extension Accessory: 83"
(211 cm)

Longitudinal Float 16" (41 cm)

Lateral Float 8" (20 cm)

Table Top Aluminum Equivalency 0.6 mm (Typical)

Table Top Length 92" (234 cm)

Table Top Width 21" (53 cm)

Trendelenber/Reverse Trendelenberg 16°

Lateral Tilt (left and right) 16°

Table Height Range Minimum 32" (81 cm)

Maximum 42" (107 cm)

Powered Floor Locks Standard

Table Shipping Weight 990 lbs (449 kg)

Patient Weight Capacity 500 lbs (226 kg)