



The **Stille Sonesta 6210** Video Fluoroscopy Table is designed for optimal imaging and translucency for procedures and examinations where a C-arm is used. The table has been developed to ensure optimal patient comfort while focusing on facilitating the positioning of the patient safely and precisely during video fluoroscopy and/or urology examinations and procedures.

Features

- Versatile – 8 positioning functions
- Maneuvered with 3 motors (height, back and seat position)
- Standing position (adjusted by back motor, maneuvered using the hand control)
- Rotation of seat/back with locked position (adjusted by back motor, maneuvered with hand control).
- Accessible for C-arm examinations
- Four double wheel rolling casters offer unencumbered maneuverability



Specifications

Dimensions

Length without foot rest/leg extension: 63 In (156 cm)

Length with foot rest/leg extension: 78 In (196 cm)

Width (total): 39 In (99 cm)

Height with seat cushion (min.): 20 In (50 cm)

Height with seat cushion (max.): 51 In (130 cm)

Imaging area without foot rest/leg extension (length): 63 In (156 cm)

Imaging area with foot rest/leg extension (length): 78 In (196 cm)

Weight capacity: 551 lbs (250 kg)

Power

120/230v 24v DC Output

Meets EMC (Electromagnetic Compatibility) for IEC 60601

All products carry the CE mark

Common Procedures

Urodynamics

Pelvic Muscle Rehabilitation

Anorectal Manometry

Pelvic Exam

Cystoscopy

Ultrasound

Vasectomy

Colposcopy

Biopsy

Colonoscopy

Gynecologic Procedure

