



Storz is one of the pioneers in the field of endoscopic imaging. Enhanced sensor technology with higher resolution and innovations in display technology as well as high-performance image processing electronics leads to continuous further developments. The **Storz Image 1 Multi-Specialty Video Tower** includes; SCB Camera Control Unit with S3 Universal Head and Coupler, Storz Xenon Light-source, 24-Inch Wide View Flat Panel Monitor, AIDA DVD Image Management System, and Storz Universal Tower Cart with LCD Monitor holders.

Components

- SCB Camera Control Unit
Comes with S3 Universal Head and Coupler.
- Storz Xenon 20133020 Light Source
- 24-Inch Wide Flat Panel Monitor
- AIDA DVD Image Management System
- Storz Universal Tower Cart
Comes with an LCD Monitor Holder



Features

Flat Panel Monitor

- 24-Inch Wide View



Storz SCB Image 1 Hub 222010 20

- Full HD Imaging System
- 1,920 x 1,080 Resolution
- Lag-Free Imaging
- Natural Color Reproduction
- Storz Image 1 S3 Universal Camera Head and Coupler



Storz Xenon 300 20133020 Lightsource

- 300 watt xenon bulb
- Video in/out BNC sockets
- 6-pin Mini-Din-Sockets for optional SCB interface



Storz AIDA DVD-M Image Management System

- Digital storage of still images, video sequences, and audio files
- Archive on DVD, CD-ROM, USB Stick, or network path
- SDI, S-Video, and composite video inputs



Storz GoKart 9601 Universal Tower Cart

- LCD Monitor Holder
- Adjustable Shelves
- Multiple Configurations

