




The standard FULL HD camera control unit with both HD and SD head receptacle. Augmented with ICM module, the unit allows printing as well as still image, movie capture and out of the box image printing by a simple push of a camera head button to suitable USB devices. The CCU is ready for OR integration through its SCB interface. Optionally the unit can be equipped with an SD-SDI interface for long-range video signal transmission.

## Features

- 1080p Full HD resolution for razor-sharp pictures without any latency
- Natural colour reproduction for excellent image rendition
- Small and compact camera head concept
- Programmable head buttons for image capture, white balance, etc.
- 3 chip Full-HD CCD technology built in
- Resolution:-1920\*1080 pixels 16:9
- Progressive scan display:-50/60 full frame per second
- Colour system:-pal/ntsc



# Specifications

## **Dimensions**

Height: 3.5 in (8.9 cm)  
Width: 12 in (30.5 cm)  
Depth: 12.5 in (31.5 cm)  
Weight: 6.8 lbs (3.1 kg)

## **Power**

Fuse: 2 x T1.6A 250V  
Power Supply Requirements: 100 - 240 V  
Power Frequency: 50/60 Hz

## **Operating and Storage Conditions**

Operating Temperature: 41° - 95 °F (5° - 35°C)  
Storage ambient Temperature: 14° to 140° F (-10° to 60°C),/li>  
Storage Humidity: 20 to 80% (non-condensing)  
Storage Atmospheric Pressure: 50 kPa to 106 kPa

