



The Karl Storz 20133120 Cermox xenon lamp by Excelitas is the xenon light source of choice for fiber optic illumination applications such as endoscopy and surgical headlamps. This is a 300 watt xenon arc surgical endoscope medical bulb with a parabolic reflector. The PE300BF light bulb has a 1000 hour life expectancy and a warranty period for up to 500 hours.

## Features

- 300 Watt Xenon bulb
- video in/out BNC sockets
- 6-pin Mini-Din-sockets for optional SCB interface



# Specifications

## **General**

Watts: 300W

Volts: 14V

Amps: 21.43A

Length: 1.97 inches - 50.0 Millimeters

