



The **Stryker 1068 Stretcher** offers the greatest head-end clearance in the industry for superior surgeon access with ample leg clearance to help facilitate surgical precision. It was designed to support delicate surgical procedures with optimal accessibility for surgeon and comfort for the patient.

Features

- Four I.V. receptacles
- Four-wheel steel ring central locking brakes with foot-end brake & steer pedal
- Dual side pedestal hydraulic controls on sides
- Hands-free Trendelenburg/Reverse Trendelenburg



Specifications

Dimensions

Overall Length: 88 In (223.5 cm)
Overall Width: 31.5 In (50 cm)
Weight Capacity: 500 lbs (228 kg)
Height Range (To Litter Top)
High: 36 In (91.4 cm)
Low: 22.5 In (57.2 cm)

Fowler Angle

0° to 90°

Knee Gatch Angle

0° to 35°

Trendelenburg / Reverse Trendelenburg

-15° to 35°

More Features

Fold-down collapsible side rails
Integrated storage tray with oxygen bottle holder
Rubber bumpers at foot-end
Stationary foot-end litter
Manual crank control backrest

