



The **Stryker 1069 Stretcher** is Ideal for either superior or temporal approach, positioning controls tuck neatly within the headpiece and allow full articulation. Enhanced-clearance and dual-articulating headpiece maximizes surgical access. Flexible tube doubles as air delivery and drape support. Dual cylinder pneumatic backrest, Dual pedestal hydraulics and Dual side hydraulic controls. The 1069 Eye Surgery Stretcher also features Fold-down Siderails, Four IV receptacles, and Four-wheel, steel-ring brakes with foot-end brake & steer pedal.

## Features

- 3" Enhanced Comfort Mattress
- 8" Omni Surface™ Casters with Wheel Covers
- Dual Cylinder Pneumatic Backrest
- dual Pedestal Hydraulics
- Dual Side Hydraulic Controls
- Fold-down Siderails
- Four IV Receptacles



## Specifications

### Dimensions

**Overall Length:** 90 In (229 cm)  
**Overall Width (Siderails Up):** 31.5 In (50 cm)  
**Weight Capacity:** 500 lbs (228 kg)  
**Height Range (To Litter Top)**  
High: 35.5 In (90 cm)  
Low: 22 In (56 cm)  
**Patient Surface:** 26 In x 87 In (66 cm x 221 cm)  
**Siderails:** 13 In x 55 In (33 cm x 140 cm)  
**Caster Diameter:** 8 In (20 cm)

### Litter Positioning

**Backrest:** 0° - 90°  
**Knee Gatch:** 0° - 30°  
**Trend. / Reverse Trend.:** ±18°

### Headpiece

**Height Adjustment:** ±4.5" (11 cm)  
**Headpiece Tilt:** ±45°

### More Features

Four-wheel, steel-ring brakes with foot-end brake & steer pedal  
Hands-free Trendelenburg  
Integrated storage tray with oxygen bottle holder  
Retractable fifth wheel  
Rubber bumpers  
Stationary foot-end litter  
Unilower pedal

