



Stryker's® eye surgery stretcher offers the greatest head-end clearance in the industry, giving surgeons ample legroom to maneuver while performing delicate eye procedures. The reliable work surface and mobile design combine to increase patient security and surgical throughput. It offers head-end clearance for surgeon access with positioning controls tucked within the headpiece to allow for full articulation.

Features

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface™ casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls Fold-down side rails
- Four I.V. receptacles
- Central-locking brakes with foot-end brake & steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece



Specifications

Length x Width 90" (229 cm) 31.5" (80 cm)

Weight Capacity 500 lbs. (227 kg)

Height Range 34" (86.4 cm)-22.25" (56.5 cm)

Backrest 0° - 90°

Trend./Reverse Trend ±18°

Siderails 13" x 55" (33 x 140 cm)

Patient Surface 26" x 87" (66 x 221 cm)

