



The **Stryker 1089 Eye Stretcher** was designed to provide comfort for the patient while offers ample accessibility for the surgeon. With its dual articulating headpiece makes it easy put the head section in any position that may be required for the procedure. The 1089 Eye stretcher 's work surface firmly locks and it's Mobile design helps to increase the patient security.

Features

- Enhanced head-end clearance for supreme surgeon access.
- Ideal for either superior or temporal approach.
- Positioning controls tuck neatly within the headpiece and allow for full articulation.
- Optional flat or contoured headpiece to comfortably support the patient during procedures.
- 500-pound weight capacity.



Specifications

Dimensions

Height: 22.25 - 39" (56.5 - 86.4 cm)

Width: 31.5" (80 cm)

Length: 90" (229 cm)

Weight Capacity: 500 lbs (227 kg)

Movements

Height Range: High: 39" (86.4 cm); Low: 22.25" (56.5 cm)

Backrest: 0° - 90°

Knee Gatch: 0° - 30°

Trend/Reverse Trend: ±18°

Headpiece Height Adjustment: ±4.5" (11 cm)

Headpiece Tilt Adjustment: ±45°

Patient Surface

Length: 87" (221 cm)

Width: 26" (66 cm)

Mattress: 3" (8 cm) Enhanced comfort mattress

Standard Features

8" (20 cm) Omni Surface casters with wheel covers

Dual cylinder pneumatic backrest

Dual pedestal hydraulics

Dual side hydraulic controls

Fold-down siderails

For IV receptacles

Central-locking brakes

Hands-free Trend.

Integrated storage tray with oxygen bottle holder

Retractable fifth-wheel steering

Rubber bumpers

Stationary foot-end litter

Uni-lower pedal

Dual-articulation, enhanced-clearance headpiece

