



The **Stryker 5050** is a multifaceted stretcher that offers a wide variety of positions. This compact stretcher chair keeps the patient's safety in mind - with hydraulic jacks and a central-locking brake for stability. The Stryker 5050 can easily be positioned from completely horizontal to a chair position to maximize efficiency and patient comfort.

Features

- Dual hydraulic jacks optimize stability for greater comfort and a more secure treatment platform.
- Pneumatic backrest that can maneuver from 0 to 90 degrees
- Dual Hydraulic jacks and central locking breaks
- 24 inch wide patient surface with a lumbar support
- Swing-down side rails with arm rests
- Hands free Trendelenburg



Specifications

Overall Length	76" (193 cm).
Width	Overall: 30.0" (76 cm). Litter: 27.5" (70 cm).
Weight Capacity	400 lbs (182 kg).
Height Range	High: 33.5" (85 cm). Low: 22.0" (56 cm).
Litter Positioning	Backrest: 0 - 90°. Leg: 0 - 80°. Trendelenburg/Reverse Trendelenburg: $\pm 45^\circ$.
Siderails	$\pm 45^\circ$. Standard: 10" x 31" (25 cm x 79 cm).
Patient Surface	24" x 74" (61 cm x 188 cm).
Caster Diameter	6" (15 cm).