



The **Stryker 940** cast cutter is a very quiet cast saw with a variable 2 speed control. Standard speed is powerful enough to remove most standard casts. High speed uses and increased speed to remove thick body jackets or hip spicas. It can be used directly plugged into a wall outlet, or connected to a vacuum like the Stryker 986 CastVac. This Stryker Cast Cutter can be used with ion nitrided blades, stainless steel blades, or titanium nitrided blades.

Features

- Variable 2 speed control.
- Can be used with vacuum or plugged into the wall.
- Quick connect vacuum mount.
- Ergonomic design.

