



The **Stryker 986 CastVac** is mounted on a five wheel mobile stand to provide easy mobility and stability. It's very quiet compared to other cast vacuum systems. The detachable canister lid makes it easy to empty the dust collection filter. The Stryker CastVac can be used with the Stryker 940, Stryker 840, and Stryker 848 cast cutters. This cast removal system features an easy-to-operate 3-way power switch.

## Features

- Removable canister lid.
- Swivel hose mount.
- Handpiece control. Vacuum automatically turns on when cast cutter, or cast saw is activated.
- Quiet vacuum.
- 2.4m vac hose.
- Dimensions: 20.1" x 20.1" x 34.6" (51 x 51 x 88 cm).

