

Features

-
- A manual height-adjustable table with a blue frame and a black padded top. The table is shown in its fully extended position. The frame is labeled 'STUOL 114' and 'STUOL 114'. The table is mounted on a blue base.

Specifications

Max. Patient weight	454 kg for normal orientation 227 kg for reverse orientation
Height adjustable (without pad)	635–1075 mm
Longitudinal slide	430 mm 1, 330 mm 2
trendelenburg /reverse Trendelenburg	30°/31°
Tilt Left/ right	+20°/–20°
Back plate up/down	+90°/–45°
Leg plate up/down	+30°/–105°
Split leg plates up/down	+30°/–90°
Head rest up/down	+45°/–90°
Kidney elevator (optional)	Lift to 75 mm
„Flex“ (Auto position)	225°
Table length	2124 mm
Tabletop width (inc. side rails)	534 mm (587 mm)
Image width	465 mm
Manual backup system	integrated
Collision protection with all powered tabletop functions	integrated
Battery capacity	5-8 working days
Segments of OR tabletop	up to 8
Tabletop material	Carbon fibre/pheno
Base material	Reinforced fibre plastics/ABS polymer, stainless steel