



Stryker's GoBed II Med/Surg bed was developed with extensive input and feedback from nursing professionals to add efficiency and mobility to patient care practices. The GoBed II is a motorized, fully-electric equipped bed with numerous one-button or one-handed controls allowing caregivers to spend their time helping patients recover quickly in a secure and comforting environment. This bed features full length siderail coverage to aid in patient safety, a low bed height of 14.5 inches promoting proper patient positioning for safe entry and exit, integrated pump holder conveniently located on footboard, and 120-volt outlet to help maximize space at the foot end of the bed. This bed is ideal for healthcare facilities for Med/Surg purposes, educational and training facilities and residential settings.

Features

- 14.5-inch (37 cm) low-bed height
- Retractable fifth-wheel steering
- Trend. display on footboard (only available with scale option)
- Centrally-located steer and four-wheel brake mechanism
- StayPut® retractable bed frame design
- 6-inch (15 cm) casters
- Trend./reverse Trend.
- Four independent electric motors
- Four drainage bag hooks
- Eight IV pole/traction equipment sockets
- Roller bumpers
- Manual back-up for head and knee control
- Integrated pump holder
- Nurse controls on footboard and siderails
- Fixed patient controls on siderails
- Degree indicator for head elevation



Specifications

Overall length	94.25" (239.4 cm)
Overall Width	Side Rails Up: 40" (101.6 cm). Side Rails Down: 39" (99.1 cm).
Weight Capacity	500 lbs (228 kg).
Height Range (to litter top)	High: 29" (73.7 cm) Low: 14.5" (36.8 cm)
Litter Positioning	Backrest: 0°– 60° (siderails down) Knee Gatch: 0°– 28° Trendelenburg/Reverse Trendelenburg: ±14° Retraction: 10" (25 cm)
Patient Surface	36" x 84" (91 x 213 cm)
Electronics	Volts: 100 VAC, 120 VAC, 200 VAC, 220 VAC Current Leakage: < 100 microamperes Frequency: 50 – 60 Hz Hospital Grade Plug: Yes.