



The Stryker MX-PRO R3 is light, durable, and fast-handling. This manual cot features a 650 lb weight capacity, yet it weighs only 81 lb. You get a lightweight cot that's strong, durable, and easy to handle. One-hand operation of backrest, side rail, breakaway head section and safety bar frees the user to perform life-saving tasks. Foot-end lift and grip design and duplicate foot-end controls provide ergonomics that accommodate different-sized operators.

## Features

- Scientifically optimized lift bar design
- Positive action height adjustment
- One-hand release breakaway head section with safety bar
- Oversize wheels maximize mobility on all types of terrain
- Wheel lock option helps keep unoccupied cot from rolling
- Sealed fluid resistant bolster mattress
- One-four point shoulder harness restraint
- Two lap belt restraints
- One-hand release, fold down side rails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- 650 pound patient capacity
- Fold-down siderails enable caregiver patient access even in the smallest spaces.
- X-frame that features 7 height positions
- Floor-mount safety hook



# Specifications

**Height**

Position 1: 13.5" (34cm)  
Position 2: 21.0" (53cm)  
Position 3: 25.5" (65cm)  
Position 4: 29.0" (74cm)  
Position 5: 32.0" (81cm)  
Position 6: 35.0" (89cm)  
Position 7: 37.5" (95cm)

**Length**

Standard: 80.5" (205 cm)  
Minumum: 62.0" (158 cm)

**Width**

23"(58 cm)

**Weight**

81 lbs (38kg)

**Wheels**

Diameter: 6" (15 cm)  
Width: 2" (5cm)

**Articulaton**

backrest: 2 - 73°  
Shock Position: + 14°

**Maximum Weight Capacity**

650 lbs (295 kg)

**Minimum Operations  
Required**

Occupied Cot: 2  
Unoccupied Cot: 1

**Recommended Loading  
Height**

Up to 32" (81cm)