



The **Stryker Power-PRO XT** is an innovative battery-powered hydraulic system that raises and lowers the patient with the touch of a button. The industry-leading Stryker Power-PRO XT powered ambulance cot dramatically reduces strenuous lifting and the associated risk of back injury. Medics experience frequent spinal loading due to repetitive motions such as lifting, lowering, carrying, and bending. Use of the Power-PRO XT has proven to reduce spinal loading, resulting in reduced injuries, lost or modified workdays, and workers' compensation costs, and increased recruitment and retention. The innovative battery powered hydraulic system raises and lowers patients with the touch of a button. The Stryker Power-PRO XT sets the standard for superior performance and outstanding ergonomics. Easy to use manual back-up system provides complete cot operation in the event of power loss. Foot end mounted Stryker SMRT battery is accessible and easily changed in any situation. Retractable head section allows cot to be shortened in any height position for maximum versatility. Exclusive automatic high speed retract feature reduces load and unload times. The Stryker Power-PRO XT – now the industry standard of care.

Features

- Innovative battery-powered hydraulic system raises and lowers the cot at the touch of a button.
- Easy-to-use manual back-up system.
- Retractable head section enables navigation in tight spaces in any height position.
- Automatic high-speed retract feature.
- XPS Siderails.
- Steer-Lock.
- SMRT Battery Kit.



Specifications

Weight Capacity	700 lbs (318 kg).
Maximum Unassisted Lift Capacity	500 lbs (225 kg).
Backrest Articulation/ Shock Position	0 to 73° / + 15°.
Overall Length	81" (206 cm).
Minimum Length	63" (160 cm).
Width	backrest: 2 - 73° Shock Position: +14°
Height	650 lbs (295 kg)
Weight	125 lbs (57 kg).
Caster Diameter	6" (15 cm).
Caster Width	2" (5 cm).
Minimum Operators Required for Loading/ Unloading	Occupied Cot: 2. Unoccupied Cot: 1.
Recommended Loading Height	Up to 36" (91 cm).
Recommended Working Height	(excluding mattress): 15.75" (40 cm).

