



The **Stryker Prime Series** is a user friendly stretcher. It was designed to reduce the strain the hospital worker while moving and still provide substantial patient comfort at the same time. The Stryker Prime stretcher offers advanced mobility with larger wheels and a zoom motorized drive system. The stretcher comes with electrical functionality making it easier for the patient to make adjustments to improve comfort without calling for a caregiver. With an electric power lift, you can simply touch the hydraulic pedal and raise the litter. The electrical system attached to the Stryker Prime can reduce the risk of the patient falling with two alarm zones. Using the system you can even weigh the patient with an integrated scale in both pounds and kilograms. The stretcher can hold a maximum weight of 700 lbs.

## Features

- Four wheel, center-locking brake system.
- Dual side-mounted foot control hydraulics with uni-lower pedal.
- $\pm$  Trendelenburg.
- Less than 21-inch low height.
- Power-washable design.
- 4 integrated I.V. receptacles.
- Glidaway siderails



## Specifications

### Dimensions

Overall Length:  
1105, 1115: 85 in. (216 cm)  
1125: 86 in. (218 cm)  
Overall Width:  
Siderails Up: 34 in. (86.4 cm)  
Siderails down: 30.25 in. (76.8 cm)  
Weight Capacity: 700 lbs. (318 kg)  
Height Range - High: 34 in. (86.4 cm)  
Height Range - Low (1105, 1115): 20.75 in. (53 cm)  
Height Range - Low (1125): 23 in. (58 cm)

### 26" and 30" Litter Specifications

Litter Positioning  
Backrest: 0° - 90°  
Trend. Reverse Trend. ( 1105, 1115): +/- 17°  
Trend. Revers Trend. (1125): +/- 16°  
Knee Gatch: 0° - 40°  
Maximum Angle: 87.5°  
Patient Surface: 26 in. x 75.25 in. (66 x 191 cm)  
Side Rails: 14 in x 57.75 in. (36 x 147 cm)  
Caster Diameter: 8 in. (20 cm)

