



The **Stryker SDP1000 Printer** is designed specifically for applications where image quality, durability, and accurate color reproduction are essential. With the combination of a 301 dpi thermal head and a built-in protective lamination feature, the **SDP1000** creates durable, printed images in full HD resolution. The SDP1000 is compatible with Stryker Digital Capture devices and the compact design and front loading features of the printer make it ideal for use in Stryker hospital carts.

Features

- Accurate and consistent color reproduction with dye sublimation technology.
- Self-laminating roll media—resistant to light, heat, water, and fingerprints.
- User-friendly front loading tray.
- Low Paper/Ink Indicator.
- Designed for compatibility with current Stryker Digital Capture Devices.
- Compact, medical grade color printer.
- Easy to operate.

