



The Stryker Secure II bed frame is a full featured and flexible bed frame that is fully customizable. The Secure II offers many features that provide the patient and caregiver with ease of care and comfort.

Responsive engineering over time has resulted in a med/surg bed that offers the highest level of security, intuitive operation and flexible functionality. Fully customizable through a host of features designed around patient and caregiver satisfaction, the Secure II represents the benchmark in patient handling.

## Features

- Dual Pedestal design.
- CPR Release.
- In-rail Fowler/Knee Gatch control.
- Electronic function lock-out control.
- Four-wheel, steel-ring brake system with centrally-located activation pedals.
- Backrest angle indicator.
- Locking-caster steering with four traction sockets.
- Removable headboard with CPR capability, retractable litter top, and four foley bag hooks.
- Four IV receptacles.
- Patient restraint locations.
- Trendelenburg/Reverse Trendelenburg controls.



## Specifications

<b>Overall length</b>	93" (236 cm)
<b>Overall Width</b>	Side Rails Up: 42.5" (108 cm). Side Rails Down: 40" (102 cm).
<b>Weight Capacity</b>	500 lbs (228 kg).
<b>Height Range (to litter top)</b>	High: 30" (76 cm). Low: 16" (41 cm).
<b>Litter Positioning</b>	Backrest: 0 - 60°. Knee Gatch: 0 - 40°. Trendelenburg/Reverse Trendelenburg: ± 12°. Retraction: 14".
<b>Patient Surface</b>	35" x 84" (89 cm x 213 cm).
<b>Electronics</b>	Volts: 115 VAC (Domestic) / 230 VAC (International). Ampere Rating: 7A (Domestic) / 4A (International) / 10A (rating for 115 volt outlet option). Current Leakage: 100 microamperes. Frequency: 60 Hz (Domestic) / 50-60 Hz (International). Hospital Grade Plug: Yes.
<b>Caster Diameter</b>	6" (15 cm).