



The Trio is a fusion of surgical stretcher, mobile operating table, and recovery platform in one unit. With the Stryker Trio the patient remains on the same surface for pre-op, surgery, and post-op. By doing so this greatly reduces patient handling, and improves patient and staff safety. All while increasing case efficiency turnover time.

Features

- 3" (8 cm) Ultra Comfort mattress with sealed seams
- 6" (15 cm) Omni Surface™ casters
- 78" (198 cm) patient surface
- Quadra Channel™ stabilizing lift system
- Siderails enable fluoro clearance
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activators
- Backrest raises to 80°
- Removable foot section (0°- 50°)
- 20° Trend./15° reverse Trend.
- Left and right lateral tilt $\pm 15^\circ$
- Ultraview litter top for unobstructed fluoroscopic access
- Three-sided foot controls for lift/lower functions
- Retractable fifth-wheel steering
- Low profile base with thermoformed hood
- Full-length surgery accessory rails



Specifications

Dimensions

Overall Length: 80" (203 cm)

Overall Width:

Siderails Up: 33" (84 cm)

Siderails Down Or Without: 30" (76 cm)

Weight Capacity: 500 lbs (228 kg)

Height Range: (to litter top)

High: 40" (102 cm)

Low: 26" (66 cm)

Litter Positioning:

Backrest: 1°- 80°

Foot Section: 0°-50°

Trend./Reverse Trend: +20°/-15° Trend

Lateral Tilt: ±15°

Patient Surface: 27" x 78" (69 x 198 cm)

Siderails: 12" x 49" (30 x 124 cm)

Caster Diameter: 6" (15 cm)

