

Features

- g
with dual side and head
- 
- The image shows a Gynnie mobile gynecological examination table. It features a white base with a control panel on the left side, including a power switch and a hand control. The table has a purple cushioned surface with a headrest and a footrest. The base is equipped with four casters for mobility. A vertical stand with a red emergency stop button is located at the head of the table.

Specifications

Dimensions

Overall Length: 81.5" (207 cm)

Overall Width

Siderails Up: 31.5" (80 cm)

Siderails Down: 26.75" (68 cm)

Weight Capacity: 500 lb. (227 kg)

Patient Surface: 26" x 77" (66 x 196 cm)

Siderails: 13" x 39" (34 x 99 cm)

Caster Diameter: 8" (20 cm)

Break away point from wall: 64" (162.6 cm)

Height Range and Positioning

Height Range (to litter top)

High: 35" (88.9 cm)

Low: 22" (55.9 cm)

Litter Positioning

Backrest: 0°– 90°

Knee Gatch: 0°– 90°

Trend./Reverse Trend.: ±18°

Maximum Angle: 87.5°



SOMA TECHNOLOGY, INC.™
New, Demo & Refurbished Medical Equipment