



The **Stryker 1710 Renaissance Stretcher** is available in four different models. Each stretcher comes with a durable patient handling system. The Standard fifth wheel steering and pop-up handles allow for easy maneuverability. Steel-ring brakes with optional four-sided brake controls ensure a firmly locked working surface.

Features

- 3" (8 cm) Ultra Comfort mattress with sealed seams
- 6" (15 cm) Omni Surface™ casters
- 78" (198 cm) patient surface
- Quadra Channel™ stabilizing lift system
- Siderails enable fluoro clearance
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activators
- Backrest raises to 80°
- Removable foot section (0°- 50°)
- 20° Trend./15° reverse Trend.
- Ultraview litter top for unobstructed fluoroscopic access
- Three-sided foot controls for lift/lower functions
- Retractable fifth-wheel steering
- Low profile base with thermoformed hood
- Full-length surgery accessory rails



Specifications

Dimensions

Overall Length: 84" (213 cm)

Overall Width:

Siderails Up: 34.5" (88 cm)

Siderails Down: 34.5" (88 cm)

Weight Capacity: 500 lbs (228 kg)

Height Range: (to litter top)

High: 35.5" (90 cm)

Low: 22" (56 cm)

Litter Positioning:

Backrest: 0°- 90°

Knee Gatch: 0°-30°

Trend./Reverse Trend:t: ±18°

Patient Surface: 29" x 76" (74 x 193 cm)

Siderails: 13" x 55" (33 x 140 cm)

Caster Diameter: 8" (20 cm)

