



A cost effective basic fixed-height Stryker 720 Transport Stretcher with Fowler (manual), Side Rails, IV Pole and O2 Tank Well.

The simple design of the Stryker Transport Stretcher provides safe and efficient transport. Ideally suited for the ASC environment. The Stryker 720 also has Central Brake and Steering functions, that features a fifth wheel mobility for easier traction and cornering. Additionally, the Stryker Transport Stretcher also features a Modern Design and a 500 Lbs Weight Capacity.

Features

- Fixed Height
- Manual Trendelenburg
- 3" (8 cm) mattress
- 8" (20 cm) Omni Surface™ casters
- 24" (61 cm) wide patient surface
- Fold-to-foot siderails
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual end activators
- Integrated storage tray with oxygen bottle holder
- Locking caster steering
- Roller bumpers
- Stationary foot-end litter
- Dual pedestal design



Specifications

Dimensions	Overall Length: 83in (211 cm) Overall Width: 30 in (76 cm) Weight Capacity: 500 lb (228 Kg)
Fixed height	28 in, 30 in or 32 in (71, 76, 81 cm)
Litter Positioning Backrest	0° - 90°
Litter Positioning Trend./Reverse Trend.	12° Trend.
Patient Surface Standard	24 in x 76 in (61 cm x 193 cm)
Patient Surface Optional	26 in x 76 in (66 cm x 193 cm)
Siderails	13 in x 55 in (33 cm x 140 cm)
Transfer Board Long	10 in x 64.5 in (25 cm x 164 cm)
Transfer Board Short	8 in x 44 in (20 cm x 112 cm)