

Stryker 782 Camera System



Features:

- Natural colors
- Highest resolution for the clearest picture possible
- Imaging System: 768 x 494 pixels
- Scanning System: Progressive Scan of 525 lines, 2:1 interlace
- Resolution: Horizontal: 600 lines, Vertical: 450 lines
- Internal/External synchronization (automatic selection)



Typical
Manufacturer's
Image

Soma Technology Inc. • 166 Highland Park Drive • Bloomfield, CT 06002 • USA
Phone: 1.800.GET.SOMA • www.SomaTechnology.com • Email: soma@somatechnology.com

