



The New TITAN “RCS” Rear Cranial Support Headband provides enhanced comfort due to the light weight and flexible material. The new front to back crown support design allows for better module adjustment and stabilization. These features help minimize the fatigue that sets in during lengthy surgical procedures.

## Features

- Lightweight, Comfortable and the new “RCS” helps minimize fatigue.
- BRIGHTEST Illumination Available.
- Headlight cables available in bifurcated and straight configuration.
- Fused fibers on the proximal end increase transmission efficiency by 10%.
- UltraGrip Handle eases cable removal from the lightsource and aids against fiber damage.
- Newly Designed Lightweight Headband with Rear Cranial Support for added comfort.
- Headlight Module with 20mm-110mm Variable Spot @ 16” Also available in 20mm-140mm Variable Spot.
- Fused 3.5mm Bifurcated UltraGrip Cable 9’ Long with Wolf Lensed Proximal End to optimize light output Also available in 20mm-140mm Variable Spot @ 16”.
- New Molded Junctions add Flexibility and the Reinforced Sheathing Protects the Fiber Bundle.
- Headband with Module Weight – 4.5 oz.
- Headband, Module and Bifurcated Cable Weight – 6.5 oz.

