



The **Terumo Sarns 11160** Dual Heater Cooler is a reliable, efficient patient cooling/heating system. Refurbished units available. Please contact us for pricing information.

Features

- Simple controls and high reliability
- Large capacity heater with four temperature settings
- Independent water supply and recirculating system



Specifications

Dimensions

Length 24.5" (62.2 cm)
Width 15.5" (39.4 cm)
Height 27.5" (69.9 cm)
Weight, empty 75 lbs (34.0 kg)
Weight, full 125 lbs (57.0 kg)
Water Tank capacity 8.5 gal (32 L)

Connections

Quick-connect fittings. Attachment fittings and 1/2" ID vinyl tubing are included in the accessories kit.

Performance

Maximum flow rate 6.5 gal/min (25 L/min) at zero head
Cooling mode 4 23 L/min ± 1 L/min
Cooling mode 3 7.5 - 8.5 L/min
Cooling mode 2 5.5 - 6.5 /min
Cooling mode 1 3.5 - 4.5 L/min