



The **Belmont RI-2** rapid infuser rapidly delivers warmed blood and fluid to combat hypothermia and blood loss for organ transplants and other surgical procedures. Weighing only 28 pounds, the Belmont RI-2 mounts on an IV pole and can infuse at rates between 2.5 and 1000 ml/min. It features 2 air detectors and 2 temperature probes.

Features

- Primes in about 13 seconds.
- Peristaltic roller pumping mechanism.
- Battery backup in case loss of power.
- LCD waterproof touch screen display.



Specifications

Dimensions
Height: 13.5" (34.29 cm)
Width: 7.5" (19.05 cm)
Depth: 12" (30.48 cm)
Weight: 28 lbs. (12.7 kg)

Pump
Type: Peristaltic roller
Flow Rates: 2.5-1,000 ml/min

Sensors
Air Detectors: 2
Temperature Probes (Infrared): 2
Pressure Transducer
Infusion Rate Sensor
Valve Activation Sensor
Door Open Detector

Display
Type: Touchscreen waterproof LCD
Touchscreen: 4.53" x 3.4" (11.51 cm x 8.63 cm)

Battery
Rechargeable
Automatically switches to battery upon power loss
All functions work on battery except heating

