



The **Tuttnauer EZ11 Plus** is a fully automatic Autoclave that comes with all the premium features of the EZ9 Plus, but with a larger chamber size. The color display gives the user a brighter viewing area. With the Dynamic air removal technology and closed-door drying, the Ez11 plus is more efficient. The tuttnauer's autoclaves use steam as a sterilizing agent to completely sterilize surgical and non-surgical instruments. Steam Sterilization of medical instruments is a safe and effective way to clean.

Features

- The newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing.
- Innovative EZView Multi-Color Display offers real-time cycle information that is easy to read, even from a distance. The EZView displays cycle parameters, temperature, and pressure in the chamber. During a cycle, each stage is displayed with an active color wheel and text.
- Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high-efficiency air pump circulates filtered air through the chamber for fast drying. 0.2 μ HEPA filter (replaceable).
- Comprehensive 2 year parts & labor warranty and 10 year chamber warranty.
- 11" diameter X 19.8" deep chamber and flexible loading options. The exceptionally deep 19.8" chamber works in a variety of loading options to suit any office:
 - Use the 5 full-sized trays for packs/pouches
 - Use the tray rack (included) to accommodate 4 full plus 4 half sized cassettes...no additional cassette rack to buy!



Specifications

Dimensions

Height: 17.3 in (440 mm)
Width: 20.9 in (530 mm)
Depth: 24.8 in (630 mm)
Weight: 122 lbs (56kg)

Chamber Dimensions

11in x 19.8 in (230 mm x x504 mm)
Chamber Volume 7.5 gal - 28.5 liters
Tray Dimensions
Height: 0.6 in (16 mm)
Wight: 6.7 in (169 mm)
Depth: 16.3 in (413 mm)
No. of trays: 5

Power

Current: 1400 W - 12 A
Voltage - Frequency: 120v - 50/60 Hz

