



The **Tuttnauer Elara 11** is a top-of-the-line, class B sterilizer for dental, medical and veterinary equipment. The Elara Autoclave comes with a built in vacuum pump that can eliminate air pockets from any type of load- maximizing the the efficiency of the steam penetration. The vacuum pump engages after the sterilization stage, for efficient post-vacuum drying.

Features

- Large 7.5 Gallon chamber within a small foot print. A perfect fit for standard 24 inch shelves or counters
- Provided with 5 stainless steel trays and can also hold cassettes
- Water filling from the front economizes on office space by eliminating the need to keep the top of the autoclave accessible
- Short cycle and high volume pre and post vacuum pump for perfect class B sterilization results
- 316L grade stainless steel chamber and door with electro-polish finish
- A built-in steam generator provides steam and facilitates a standby steam mode
- The door closes smoothly and locks automatically



Specifications

Dimensions

Height: 17.4 in (444 mm)
Width: 20.8 in (530 mm)
Depth: 25.4 in (645 mm)
Weight: 148lbs (67 kg)

Chamber Dimensions

11 in x 19.8 in (280mm x 504 mm)
Chamber Volume 7.5 gla - 28.5 liters
Tray Dimensions
Height: 0.8 in (22 mm)
Width: 6.7 in (170 mm)
Depth: 16.3 in (442 mm)
No. of Trays: 5

Power

Current: 2300 W - 10 A
Voltage - Frequency: 1ph 230v - 50/60Hz

